

**[GSIS]**

# **How to install iOS App**

2023. 10. 11


**SAMSUNG MEDISON**

# Contents

- I . How to install iOS App
- II. How to Log-in to App

# I . How to install iOS App - 1/4

1. Delete the existing GSIS App.  
Run the Safari or Chrome App.



2. Please access the following path through your iPhone download the installation program. ★ **New version 2.00\_6**  
<https://gsis.samsunghealthcare.com/AppDownload4.html>




3. When the installation pop-up appears, select 'Install'.



# I . How to install iOS App - 2/4

4. The GSIS 2.0 App is installed on the main screen. Run the App.



5. When the 'Untrusted Enterprise Developer' message appears, select 'Cancel'.



6. Select on 'Settings'.




# I . How to install iOS App - 3/4


7. Navigate to the path below.

[Settings > General > VPN & Device Management](#)


Select 'VPN & Device Management'.



8. Select 'SAMSUNG MEDISON '.




9. Select Trust 'SAMSUNG MEDISON '.




# I . How to install iOS App - 4/4


10. In the Trusted app pop-up, Select 'Trust'.



11. The setting is completed.




12. Run the GSIS App.




# II. How to Log-in to App


1. After running the app, enter your ID / PW and log in.



2. When you log in after running the app for the first time, you must set a phone lock password. The phone lock password is to set a password for your personal phone, and you must remember it.



3. GSIS Main page



**Thank You**